
TummyCare-Max


Simple Mixing Guide

For use with Prilosec 20mg tablets at a concentration of 3 mg / mL concentration


Your kit includes:



Step 2: Add water to bottle




Step 1: Add omeprazole to grinder



Coffee grinder such as Krups, Toastmaster or Quiseen (see below). The Quissen is rated very high and is \$9.99



Step 3: Dispense



Remember:

- Refrigerate between doses.
- Do not give with acidic foods (fruit, fruit juice, fruit yogurt).
- Use within 28 days of preparing the suspension.
- The kits are reusable. Please save the bottle, syringes, and funnel for future use. You will need only to reorder refill packets.

Recommended Dosing

Total volume	3 mg/mL
dispensed	concentration
4.0mL	12 mg
4.5mL	13.5 mg
5.0mL	15 mg
5.5mL	16.5 mg
6.0mL	18 mg
6.5mL	19.5 mg
7.0mL	21 mg
7.5mL	22.5 mg
8.0mL	24 mg
8.5mL	25.5 mg
9.0mL	27 mg

- 1mL = 1cc
- 3 milliliters (mL) is the minimum volume per dose of TummyCare Max that should be given. (This is to ensure enough buffer to neutralize stomach acid)
- Doses can be given regardless of meal-time.
- For more information on dosing, visit www.infant-acid-reflux-solutions.com.